

Storage

Store the Viromer® PLASMID kit at 4°C.

Before you start

Warm all reagents to room temperature.

Optional: Change media before transfection to remove dead cells.

Cell density and confluency

Target: 60-80% at the time of transfection

Troubleshooting

In case of toxicity

Change media 4h after transfection. Use a lower amount of transfection complex on your cells. Use less Viromer® in step 2 of the protocol.

All cells are different so a good optimization is requested each time you start working with a new cell type or plasmid. We highly recommend to do the detailed optimization experiment on page 3, 5 or 7 at the beginning.

	96-well	48-well	24-well	12-well	6-well	10cm dish
medium / well	100µl	250µl	500µl	1ml	2ml	10ml
volumes used for preparing transfection complexes	1/5 x	½ x	standard	x2	x4	x20


You will find detailed optimization guides for 96-, 24- and 6-well transfections in the following pages. If you want to work in a different well-format, up- and downscaling is done easily by having a look at the volume of medium used in the different well sizes and adapt the given protocol volumes to the larger or smaller well size.

Standard Protocol

Step 1: Preparation of pDNA
Dilute your pDNA stock solution in provided Buffer PLASMID at 11 ng/μl. Prepare a volume of 330μl.

Step 2: Preparation of Viomer®
Place a **0,2μl** droplet Viomer® onto each wall of 4 fresh tubes. Immediately add **4,8μl** of Buffer PLASMID and vortex 3-5s.
Always add buffer to Viomer® and not vice versa!


Provide Viomer® into 4 fresh tubes.



Dilute Viomer® in buffer.



Step 3: Complexation
Pipette the following volumes of pDNA (step 1) on the 5μl of Viomer® solution from step 2.



Mix swiftly and incubate 15min at room temperature.



Step 4: Add 10μl of the transfection complexes on the cells
pDNA per well: 100ng

Please perform triplicates or see next page for the detailed optimization guide!

Step 5: Read-out
Incubate cells as usual. There is no need to change medium unless high amounts of transfection complex cause toxicity.
Monitor effects 24-72 hours post-transfection and determine the best conditions for your special cells.

Like we humans, all the cell types are different. A protocol which is working great in one cell type is not always transferable to a different cell type. That's why a good optimization is necessary to achieve the efficiencies you are looking for.

This is a more detailed optimization guide where you can do the most important steps in one experiment.

- playing on the amount of pDNA on the cells
- playing on the Viomer® – pDNA ratio for a better packing

It will help you to determine the optimal conditions for your special cells and plasmid much faster.

Step 1: pDNA 11ng/μl provide 280μl

		tube				
		①	②	③	④	
Step 2:	Viomer® PLASMID in μl	0,24	0,2	0,16	0,12	Always add buffer to Viomer®, not vice versa! vortex 3-5s
	Buffer pDNA in μl	5,76	4,8	3,84	2,88	
		54	67,5	72	67,5	

Step 3: **complexation** Pipette x μl of your pDNA from step 1 into the 4 tubes with diluted Viomer® of step 2. Mix swiftly and incubate 15min at room temperature.

Transfer x μl of complexes into the wells according to the pipetting scheme.

Step 4: **Forward transfection:** Add transfection complexes onto the cells seeded a day before. Mix carefully.

Reverse Transfection: Add transfection complexes to empty wells and seed the cells (in 100μl) immediately afterwards. Mix carefully.

96-well	1	2	3	4	5	6
tube 1	2μl	5μl	10μl	15μl	20μl	
tube 2	2μl	5μl	10μl	15μl	20μl	
tube 3	2μl	5μl	10μl	15μl	20μl	
tube 4	2μl	5μl	10μl	15μl	20μl	
	20ng DNA/well	50ng DNA/well	100ng DNA/well	150ng DNA/well	200ng DNA/well	

← playing on the amount of complexes / pDNA on the cells →

↑ playing on the Viomer® - pDNA ratio ↓


Step 5: **Read-out** Incubate cells as usual. There is no need to change medium unless high amounts of transfection complex cause toxicity. Monitor effects 24-72 hours post-transfection and determine the best conditions for your special cells.

Standard Protocol


Step 1: Preparation of pDNA
Dilute your pDNA stock solution in provided Buffer PLASMID at 11 ng/ μ l. Prepare a volume of 330 μ l.

Step 2: Preparation of Viomer®
Place a **0,2 μ l** droplet Viomer® onto each wall of 4 fresh tubes. Immediately add **4,8 μ l** of Buffer PLASMID and vortex 3-5s.
Always add buffer to Viomer® and not vice versa!


Provide Viomer® into 4 fresh tubes.



Dilute Viomer® in buffer.



Step 3: Complexation
Pipette the following volumes of pDNA (step 1) on the 5 μ l of Viomer® solution from step 2.



Mix swiftly and incubate 15min at room temperature.



Step 4: Add 50 μ l of the transfection complexes on the cells
pDNA per well: 500ng

Please look also at the next page for the detailed optimization guide!

Step 5: Read-out
Incubate cells as usual. There is no need to change medium unless high amounts of transfection complex cause toxicity.
Monitor effects 24-72 hours post-transfection and determine the best conditions for your special cells.

Like we humans, all the cell types are different. A protocol which is working great in one cell type is not always transferable to a different cell type. That's why a good optimization is necessary to achieve the efficiencies you are looking for.

This is a more detailed optimization guide where you can do the most important steps in one experiment.

- playing on the amount of pDNA on the cells
- playing on the Viomer® – pDNA ratio for a better packing

It will help you to determine the optimal conditions for your special cells and plasmid much faster.

Step 1: pDNA 11ng/μl provide 1200μl

		tube				
		①	②	③	④	
Step 2:	Viomer® PLASMID in μl	1,2	0,8	0,6	0,6	Always add buffer to Viomer®, not vice versa! vortex 3-5s
	Buffer pDNA in μl	28,8	19,2	14,4	14,4	
		270	272	270	339	

Step 3: **complexation** Pipette x μl of your pDNA from step 1 into the 4 tubes with diluted Viomer® of step 2. Mix swiftly and incubate 15min at room temperature.

Transfer x μl of complexes into the wells according to the pipetting scheme.

Step 4: **Forward transfection:** Add transfection complexes onto the cells seeded a day before. Mix carefully.

Reverse Transfection: Add transfection complexes to empty wells and seed the cells (in 500μl) immediately afterwards. Mix carefully.


24-well	1	2	3	4	5	6	
tube 1	10μl	25μl	50μl	75μl	100μl		↑ playing on the Viomer® - pDNA ratio ↓
tube 2	10μl	25μl	50μl	75μl	100μl		
tube 3	10μl	25μl	50μl	75μl	100μl		
tube 4	10μl	25μl	50μl	75μl	100μl		
	100ng DNA/well	250ng DNA/well	500ng DNA/well	750ng DNA/well	1000ng DNA/well		
	← playing on the amount of complexes / pDNA on the cells →						

Step 5: **Read-out** Incubate cells as usual. There is no need to change medium unless high amounts of transfection complex cause toxicity. Monitor effects 24-72 hours post-transfection and determine the best conditions for your special cells.


Standard Protocol

- Step 1: Preparation of pDNA**
Dilute your pDNA stock solution in provided Buffer PLASMID at 11 ng/μl. Prepare a volume of 1300μl.
- Step 2: Preparation of Viomer®**
Place a **0,8μl** droplet Viomer® onto each wall of 4 fresh tubes. Immediately add **19,2μl** of Buffer PLASMID and vortex 3-5s.
Always add buffer to Viomer® and not vice versa!


Provide Viomer® into 4 fresh tubes.



Dilute Viomer® in buffer.



- Step 3: Complexation**
Pipette the following volumes of pDNA (step 1) on the 5μl of Viomer® solution from step 2.



Mix swiftly and incubate 15min at room temperature.



- Step 4: Add 200μl of the transfection complexes on the cells**
pDNA per well: 2μg
Please look also at the next page for the detailed optimization guide!

- Step 5: Read-out**
Incubate cells as usual. There is no need to change medium unless high amounts of transfection complex cause toxicity.
Monitor effects 24-72 hours post-transfection and determine the best conditions for your special cells.

Like we humans, all the cell types are different. A protocol which is working great in one cell type is not always transferable to a different cell type. That's why a good optimization is necessary to achieve the efficiencies you are looking for.

This is a more detailed optimization guide where you can do the most important steps in one experiment.

- playing on the amount of pDNA on the cells
- playing on the Viomer® – pDNA ratio for a better packing

It will help you to determine the optimal conditions for your special cells and plasmid much faster.

Step 1: pDNA 11ng/μl provide 4800μl

		tube				
		①	②	③	④	
Step 2:	Viomer® PLASMID in μl	4,8	3,2	2,4	2,4	Always add buffer to Viomer®, not vice versa! vortex 3-5s
	Buffer pDNA in μl	115,2	76,8	57,6	57,6	
		1080	1088	1080	1080	

Step 3: **complexation** Pipette x μl of your pDNA from step 1 into the 4 tubes with diluted Viomer® of step 2. Mix swiftly and incubate 15min at room temperature.

Transfer x μl of complexes into the wells according to the pipetting scheme.

Step 4: **Forward transfection:** Add transfection complexes onto the cells seeded a day before. Mix carefully.

Reverse Transfection: Add transfection complexes to empty wells and seed the cells (in 2ml) immediately afterwards. Mix carefully.

6-well	1	2	3	4	5	6
tube 1	40μl	100μl	200μl	300μl	400μl	
tube 2	40μl	100μl	200μl	300μl	400μl	
tube 3	40μl	100μl	200μl	300μl	400μl	
tube 4	40μl	100μl	200μl	300μl	400μl	
	400ng DNA/well	1μg DNA/well	2μg DNA/well	3μg DNA/well	4μg DNA/well	

← playing on the amount of complexes / pDNA on the cells →

↑ playing on the Viomer® - pDNA ratio ↓

Step 5: **Read-out** Incubate cells as usual. There is no need to change medium unless high amounts of transfection complex cause toxicity. Monitor effects 24-72 hours post-transfection and determine the best conditions for your special cells.