



[VB-01LB-00, VB-01LB-01](#)



[VG-01LB-00, VG-01LB-01](#)

Fig. 1: Knockdown of gene X expression in NCH 149 glioblastoma cells transfected with Viomer® BLUE or Viomer® GREEN

- 24 h incubation
- final siRNA concentration: **10 nM**
- Viomer concentration: 10x lower than in the standard protocol

“The results are more than satisfying, considering that we were not able to transfect these cell lines with any other transfection reagent that is on the market.”

Data from I. Dokic, DKFZ Heidelberg (Germany)